

CiD clear image devices

CR Cassette Protector



No More Broken Cassettes!

The unique design of CiD's CR Cassette Protector provides an easy way to take weight-bearing X-rays of feet (AP view) while protecting CR cassettes. The CR Cassette Protector is simply placed over a cassette which is first positioned on the floor. Patients step onto the CR Cassette Protector with as much weight as needed to get the desired image.

10 x 12 in. (23 x 30cm) models:

Item 2934: 275 lb capacity

Item 2934H: 450 lb capacity

14 x 17 in. models:

Item 3947: 275 lb capacity

Item 3947H: 450 lb capacity

Clear Image Devices' series of patient positioning devices is the result of innovative research and design, responding to the expressed needs of radiographers. CiD produces only high-quality products that enhance medical care and promote the continued advancement of medical imaging.

- *Optimal positioning solutions for weight-bearing AP foot images of patients up to 450 lbs*
- *Keeps patient's foot as close to the target image platform as possible to ensure high quality images with no significant magnification*
- *X-ray lucent faceplate allows X-rays to pass through the CR Cassette Protector with very little absorption or scattering*
- *Non-slip rubber grips keep the CR Cassette Protector securely in place on a hard, smooth floor*
- *Easy to keep clean, the CR Cassette Protector is made from germ-resistant, non-porous polyethylene, polycarbonate, and stainless steel*
- *Simple, reliable, and virtually indestructible*
- *Available for all sizes and models of film, CR, and DR imaging systems*

Clear Image Devices LLC

3930 Michael Rd. North, Ann Arbor, MI 48103

Email: sales@clearimagedevices.com

Tel.: 734.645.6459 Fax: 734.661.0220



www.clearimagedevices.com

Copyright © 2008 Clear Image Devices LLC. All rights reserved.